

# Bonding / Cementation Chart



RESTORATION TYPE	MANUFACTURE BRAND RECOMMENDED	SPECIFIC TYPE	COMMENTS
<b><i>Empress Esthetic Veneer™</i></b>	<b><i>Variolink Veneer™</i></b>	Bonded Resin	Best for esthetic type restorations. 7 try-in & value shades. Easy clean-up.
<b><i>Empress Esthetic™</i></b> Inlay/onlay & crowns	<b><i>Multilink™</i></b> also Variolink II & 3M RelyX Note: When bonding to zirconia, Multilink is <u>ONLY</u> recommended	Bonded Resin	Do <u>not</u> use Resin-onomer or any hybrid cement. Expansion is problematic.
<b><i>E-MAX™</i></b>	<b><i>Variolink II™</i></b> May also use: Relyx Unicem, Fuji 1 or Ketac Cem Note: When bonding to zirconia, Multilink is <u>ONLY</u> recommended	Bonded Resin or Resin Cement	Do <u>not</u> use Resin-onomer or any hybrid cement. Expansion is problematic.
<b><i>LAVA™ Zirconia Imagine™</i></b> or other brands	<b><i>RelyX™ Unicem Self-Adhesive Universal Resin Cement</i></b>	Resin Cement	
<b><i>Implants Metal to Metal</i></b>	<b><i>Improv™</i></b>	Acrylic Urethane Base Cement	Good for provisionals & permanent bond. Achieves maximum strength over time.
<b><i>Metal/PFM</i></b> Crown & bridge, endodontic posts.	<b><i>RelyX™ Unicem Self-Adhesive Universal Resin Cement</i></b>	Resin Cement	